



Rocketmedical

Rocket® LEEP Loops



The Rocket® LEEP Loop Range was the world's first single use loop electrode designed specially for LEEP.

20 years on the range is still in constant production. It continues Rocket Medical's tradition of providing exceptional value, reliability and quality

Designed to overcome the problems of high cost, cross infection risks and health and safety hazards commonly associated with reprocessing reusable LEEP electrodes.

- **UNIQUE 0.15mm Rocket® wire** versions to deliver higher cutting performance with significantly less thermal damage to the sample. (Small and Square styles only).
- **2.4mm (3/32") x 10cm shaft** to fit all standard UK, US and European diathermy systems
- **DELTA™** style loops for increased wire rigidity and increased control of penetration.

Rocket® Range		Size	Code
Rocket® LEEP Loop, for single patient use, packed sterile in cartons 25 units	Small	14mm	R57800
	Medium	18mm	R57801
	Large	22mm	R57802
	Square	10mm	R57803
Rocket® DELTA™ LEEP Loop, for single patient use, packed sterile in cartons 25 units	Medium	18mm	R57851
	Large	22mm	R57852
Rocket® Ball Electrodes, for single patient use, packed sterile in cartons 25 units	Standard	5mm	R57804
	Small	3mm	R57807
Rocket® Diathermy Extension, for single patient use, individually packed sterile		50mm	R57818
Rocket® NETZ™ Needles, for single patient use, packed sterile in cartons 25 units		15mm	R57860





Rocketmedical



Rocket Medical plc 150 Recreation Park Drive Unit 3 Hingham MA 02043 USA
T: 1-781-749-6223 F: 1-781-749-6235 E: usa@rocketmedical.com www.rocketmedical.com

Registered office: Rocket Medical plc Imperial Way Watford WD24 4XX England. Company Number: 03276608