


# Rocket® Marrowgauge™ Needle



Rocket® Marrowgauge™ Needles provide a highly cost effective solution to the many problems associated reprocessing of conventional marrow aspiration needles.

- **TWO PRIMARY NEEDLE LENGTHS:** to meet a wide range of patient needs.
- **CONTROLLED DEPTH** of penetration with 2 sleeves allowing sternal and iliac crest puncture with total control of penetration depth.
- **GROUND TAPER TIP** with precision ground tip to enable easy penetration of bone.
- **ERGONOMIC HANDLE** sits comfortably in the operator's hand during initial penetration.
- **SIMPLE** locking mechanism secures stylet to ensure clean controlled bone penetration.

| Ordering Information   | Sternal<br>Iliac Crest | Code          |
|--|------------------------|---------------|
| <b>Rocket® Marrowgauge™ Needles</b><br>15g, with moulded body, retaining cap, stylet and removable sleeves.<br>Peel pouch packed in cartons of 10 needle sets, sterile, for single patient use | 10-22mm<br>22-39mm     | <b>R56780</b> |
|  | 35-47mm<br>47-64mm     | <b>R56781</b> |

