



Rocketmedical

Rocket® OHSS & Cyst Aspiration Sets

Rocket® OHSS Aspiration Sets offer a fast, reliable solution for draining ovarian cysts and for the treatment of Ovarian Hyperstimulation Syndrome in the office setting.

- **HIGH RESOLUTION** 16G x 35cm echogenic needle prepared to give smoother, easier tissue penetration with clear ultrasound definition.
- **HIGH VACUUM ASPIRATION BOTTLES** designed for connection to Craft™ Suction Pumps delivering up to 400mmHg for ease of use in the office or clinic setting. Each pack contains either 2 x 250 bottles or 1 x 1000ml to minimize operating time.
- **HIGH ASPIRATION RATES:** Rocket® OHSS sets typically deliver up to 190ml/min @ 400mmHg giving rapid, controlled fluid aspiration.
- **SOFT, KINK FREE TUBING SET** dramatically improves ease of handling and eliminates the kinking problems common to rigid FEP tubing.
- **STERILE** - Individually pouch packed, intended for one time use.

Specifications:

16g x 35cm needle 15° background bevel with echogenic surface.
 Supplied with 2 x 250ml or 1 x 1000ml vacuum collection bottles and 700mm connection set.
 Suitable for use with all **Rocket® Craft™ Suction Pumps** delivering 400mmHg.

Ordering Information	Size/Length	Code
Rocket® OHSS & Cyst Aspiration Set complete with 2 x 250ml drainage bottles	16g x 35cm	R57615-16-00
Rocket® OHSS & Cyst Aspiration Set complete with 1 x 1000ml drainage bottle	16g x 35cm	R57616-16-00
Rocket® Replacement Drainage Bottle 1 x 1000ml drainage bottle		R57616-BT-00

