

# CUSTOMER INFORMATION BULLETIN

## Rocket Seldinger Kits: Change to method of catheter fixation.



### Dear Customer,

As part of our programme of continuous improvement we are delighted to announce that we are making a change to the method of catheter fixation supplied with our Rocket Seldinger Kits.

This change has been introduced following detailed customer research during recent clinical evaluations.

As a result of customer feedback, we are removing the 1.0 silk suture and needle from our packs and replacing it with a new, self adhesive drain fixation device which offers a number of important advantages:

- Reduces risk of accidental displacement of catheter.
- Safe anchoring and protection of catheter for up to one week.
- Allows inspection of puncture site.
- Reduces catheter movement and consequent enlargement of puncture site.
- Enables the patient to shower with catheter inserted.
- Prevents kinking of catheter.
- Prevents contamination of puncture site from environment.
- Available in 2 sizes: for Rocket 12 & 16Fg catheters, and for 18 & 20Fg.
- Included in our Rocket Seldinger SG & PK packs as standard.
- May be used instead of or in addition to, a suture of your choice.



Customer Information Bulletins on all of our products can be found on our web site at:

[http://www.rocketmedical.com/Rocket\\_support\\_CIB\\_main.htm](http://www.rocketmedical.com/Rocket_support_CIB_main.htm)

**If you have any questions, please contact the Customer Services Team for assistance:**

Tel: +44 191 419 6988. Fax: +44 419 6889. Email: [customerservices@rocketmedical.com](mailto:customerservices@rocketmedical.com)

