

Rocket®

Access Sets and Accessories

Rocket Medical plc offers the best in simple and multi-part introducer needles, dilators, guidewires and custom access sets.

Description		Length	Code
Introducer Needles			
IR™ Introducer Needles 5FG outer sheath, 19G inner cannula tube, 21G tri-facet trocar needle	R/O	19.5cm	R58.900
	Clear	19.5cm	R58.900-00-KL
	R/O	15cm	R58900-00-15
	Clear	15cm	R58900-KL-15
IR™ DUO™ Introducer Needle 5FG outer sheath, 19G inner chiba needle	Clear	13cm	R58900-CH-15
Seldinger Needle, sterile, 10 per pack	18G	7cm	R58.902
Chiba Needle, sterile, 10 per pack	22G	15cm	R58.903
Guidewires			
Fixed Core, PTFE coated, 3mm J tip, dispenser loaded. Sterile per 10	0.035"	80cm	R58.905-35-80
Dilators			
Radiopaque, luer lock, polyethylene, barium loaded tapered dilators, sterile per 10	4-20FG	21cm	R58.921-FG-21
Pre-packed Access Sets			
PTFE J tip wire, radiolucent dilators, 18G x 7cm seldinger needle, sterile each	0.018" 0.035" 4-8FGFG	80cm 21cm	R58.912

