

Rocket® Introducer Needles



Rocket Medical plc offer their latest developments in three and two part introducer needle sets. We have listened carefully to our users and have added unique features to meet those needs.

- **LOCKING HUB:** All **Rocket® Introducer Needles** feature a fully locking hub system to make insertion easier and more secure.
- **RADIOPAQUE & CLEAR** (Kellert style) versions are available as standard to suit all users. Outer permits up to 0.038” wire, inner cannula up to 0.032” wire.
- **SMOOTH TRANSITION** between the tri-facet or chiba needle, cannula and outer sheath to ensure trouble free insertion with the minimum of trauma.
- **SILICONISED** needle sets ensure smooth, drag free withdrawal after insertion.
- **COST EFFECTIVE:** UK designed and manufactured to minimise procedure costs.



Description	Sheath	Length	Code
Rocket® Introducer Needles 5FG outer sheath, 19G inner cannula tube, 21G tri-facet trocar needle Sterile, for single patient use, pack in carton of 10 units	R/O	19.5cm	R58900
	Clear	19.5cm	R58900-00-KL
	R/O	15cm	R58900-00-15
	Clear	15cm	R58900-KL-15
Rocket® DUO™ Introducer Needle 5FG outer sheath, 19G inner chiba needle	Clear	15cm	R58900-CH-15

