

Rocket®

Multipurpose Drainage Catheters



Rocket Medical offer a new solution to the challenge of maintaining a soft, kink free material with the need to prevent ruckle back during insertion.

- **EasyTip™ CORE** to prevent ruckle-back, minimising resistance during insertion; reducing discomfort and speeding correct placement (8.4Fg-14FG).
- **SOFT, KINK FREE** catheter material, improves patient comfort and prevents unnecessary catheter replacement caused by harder, 'kink resistant' materials.
- **COMPLETE WITH BAG CONNECTION SET** to reduce costs and make Rocket® MPD Catheters even better value for your department or ward.
- **SECURE LOCKING SYSTEM** to give reliable pigtail locking without cumbersome 'keys' or bulky hubs. Self adhesive removal instructions (for addition to the patient's notes) are provided with every locking drain.

Description	Size	Code
Rocket® Multipurpose Drainage Set: plain, pigtail catheter, trocar and stiffener set, 30cm x 6 drainage holes complete with bag connection set	6FG	R58810-06-30
	8FG	R58810-84-30
	10FG	R58810-10-30
	12FG	R58810-12-30
	14FG	R58810-14-30
Rocket® Multipurpose Drainage Set: plain coated, pigtail catheter, trocar and stiffener set, 30cm x 6 drainage holes complete with bag connection set	6FG	R58820-06-30
	8FG	R58820-84-30
	10FG	R58820-10-30
	12FG	R58820-12-30
Rocket® Multipurpose Drainage Set: locking pigtail catheter, trocar and stiffener set, 30cm x 6 drainage holes complete with bag connection set	6FG	R58830-06-30
	8FG	R58830-84-30
	10FG	R58830-10-30
	12FG	R58830-12-30
Rocket® Multipurpose Drainage Set: coated, locking catheter, trocar and stiffener set, 30cm x 6 drainage holes complete with bag connection set	14FG	R58830-14-30
	6FG	R58840-06-30
	8FG	R58840-84-30
	10FG	R58840-10-30
Rocket® Multipurpose Drainage Set: coated, locking catheter, trocar and stiffener set, 30cm x 6 drainage holes complete with bag connection set	12FG	R58840-12-30
	14FG	R58840-14-30
	5FG	R58851-05-30
	6FG	R58851-06-30
Rocket® Multipurpose Exchange Drainage Catheters, plain pigtail 30cm x 6 drainage holes	7FG	R58851-07-30
	8FG	R58851-84-30
	10FG	R58851-10-30
	12FG	R58851-12-30
	14FG	R58851-14-30
Rcket® Multipurpose Exchange Drainage Catheters, locking pigtail 30c0m x 6 drainage holes	6FG	R58856-06-30
	7FG	R58856-07-30
	8FG	R58856-84-30
	10FG	R58856-10-30
	12FG	R58856-12-30
	14FG	R58856-14-30

