


# Rocket® Multipurpose Needle Drain & Ward Procedure Tray

**Rocket® Multipurpose Needle Drain and Multipurpose Ward Procedure Trays** offer a purpose designed, cost-effective solution for routine fluid drainage procedures in the department or ward settings.

#### Indications:

For the drainage of body fluids including ascites, urine and pleural fluid. For insertion into, but not restricted to, the abdominal cavity, chest and urinary bladder using a 'direct stab' technique combined with aspiration and/or ultrasound guidance.



- **CE MARKED, PURPOSE DESIGNED DEVICE**, designed specifically for paracentesis/drainage of abdominal ascites and other drainage applications which avoids the risks of inappropriate use of non specific devices such as generic suprapubic drains.
- **COMPLETE SET**: Avoids wasted time and increased cost in searching for separate components for the procedure.
- **SOFT, KINK FREE** 8FG catheter material, improves patient comfort and prevents unnecessary catheter blockage caused by thin walled materials.
- **EasyTip™ CORE** combined with a long taper tip prevents ruckle-back and minimises resistance during insertion, so reducing discomfort and easing correct placement
- **SILICONISED NEEDLE SET** to aid needle withdrawal, improving ease of insertion and reducing patient discomfort.
- **SOFT INTEGRAL STITCH PLATE** combined with flexible tube support, provides simple, secure fixation after insertion.
- **REMOVABLE NEEDLE CAP** permits aspiration to aid correct placement during insertion.
- **DRAINAGE BAG CONNECTION SET** included to make insertion and connection to a wide range of drainage options quick and secure.

Description	Order Code	NHS Logistics
<b>Rocket® Multipurpose Needle Drain Set</b> 8FG x 20cm x 6 holes with soft stitch plate and insertion needle set. 30cm drainage bag connection set with 2 way stopcock, No.11 scalpel . Instructions for Use Sterile, for single patient use, supplied in multiples of 5 sets.	<b>R58800-08-ND</b>	<b>FSW372</b> <b>£33.99 inc VAT</b> in multiples of 5 units <b>£35.99 inc VAT</b> in singles
<b>Rocket® Multipurpose Needle Drain Procedure Tray</b> 8FG needle drain with trocar needle, 30cm drain bag connection set + 2 way tap, 2000ml drainage bag with empty port, 20ml, 12ml syringes, luer lock, 21g x 1.5" needle, 25g x 5/8" needle, 5in Sharp/Blunt Scissors, 60ml gallipot, Pair polyco gloves, medium, Gauze swabs x 5, Sponge stick, 8", Scalpel No.11, 11x15cm 'Mepore' dressing, 50x60m adhesive fenestrated drape, '0' Mersilk suture on curved needle, Instructions for Use Sterile, for single patient use, supplied individually.	<b>R58800-08-NK</b>	<b>FSW373</b> <b>£69.00 inc VAT</b> (Blue Diamond)

