

Rocket® PCA

Percutaneous Access Sets



Rocket Medical plc offer their latest development for routine and difficult percutaneous access. We have listened carefully to our customer's views and have added unique features to meet those needs. Designed and manufactured in UK, the Rocket® PCA Set offers reliable, cost effective and secure access.

- **LOCKING HUB, All Rocket® PCA Sets** feature a fully locking hub system to make insertion easier and more secure.
- **NITINOL GUIDEWIRE** 0.018" x 60cm with 8cm soft platinum alloy tip for rapid, controlled wire positioning.
- **HIGH DENSITY RADIOPAQUE** outer sheath as standard, RO banded version also available for increased visibility.
- **SMOOTH TRANSITION** between the inner dilator and outer sheath to ensure trouble free insertion with the minimum of trauma.
- **SILICONISED** needle sets ensures smooth, drag free withdrawal after insertion.
- **COST EFFECTIVE:** UK designed and manufactured to minimise procedure costs.



Description	Size	Code
Rocket® Percutaneous Access Sets 7FG outer sheath, 4FG inner dilator, 22G stiffening cannula, length: Nitinol Guidewire, 0.018" x 60cm, 22G x 15cm Chiba Needle	7FG x 19cm	R58913

