

Rocket® Bone Marrow Biospy Kit

The Rocket® Bone Marrow Biospy Kit is the direct result of recent user research. The new design includes all those items most commonly identified as required in the ideal bone marrow aspiration kit.

- **COMPREHENSIVE KITS:** designed to contain the most commonly required items for fast, reliable bone marrow aspiration and marrow biopsy.
- **REDUCES RISK AND STRESS** caused by searching for common items in the ward or department and struggling with missing or worn instruments which are frequently common in reusable trays.
- **KIT CONTENTS: Sterile, for single patient use, supplied in cartons of 10 packs**
- **REDUCES WASTE:** Combining components into a single kit significantly reduces wasted packaging and reduces wasted time resulting in manufacturing cost savings which means a better utilisation of Trust resources.



■ **KIT CONTENTS:**

- 1 x Rocket Marrowgauge Bone Marrow Aspiration Needle 15G
- 1 x Rocket Trephine Bone Marrow Biospy Needle 11G
- 3 x 30ml Universal Containers, capped with labels
- 10 x VFM microscope slides, twin frost, square edge
- 2 x 5 slot microscope slide transport boxes
- 3 x 2ml Syringe luer lock
- 1 x 12ml Syringe luer lock
- 1 x 20ml Syringe luer lock

- 1 x Microlance 21g /1.5" green needle
- 1 x Microlance 25g 5/8" yellow needle
- 5 x 10cm x 10cm gauze swabs
- 1 x 60ml gallipot
- 1 x 8" skin prep sponge stick
- 1 x fenestrated drape
- 1 x square drape 1m²
- 1 x 5in sharp/blunt scissors
- 1 x Mepore wound dressing 110mm x 150mm

Ordering Information	Code
<p>Rocket® Bone Marrow Biopsy Kits Packed in ward procedure tray over wrapped with field drape, in cartons of 5 kits sterile, for single patient use</p>	<p>R56750</p>

