

Rocket®

Indwelling Pleural Catheter

The Rocket® Indwelling Pleural Catheter or IPC range is indicated for intermittent, long term drainage of symptomatic, recurrent, pleural effusion. Including malignant pleural effusions that do not respond to medical treatment of underlying disease.

For further information please visit our website: www.rocketmedical.com or contact our Customer Services Team to arrange a demonstration from one of our Clinical Sales Specialists.

The system is indicated for:

1. The palliation of dyspnoea due to pleural effusion.
2. Providing pleuradhesion (resolution of the pleural effusion)

The system features:

UNIQUE FLOW CONTROL:

- Gives patient or carer total control over fluid flow, significantly reducing discomfort during drainage.
- Reduces risk of accidental vacuum release prior to catheter connection.

SILICONE CATHETER:


- Soft, kink resistant material approved for long use, giving increased patient comfort.
- 16Fg x 40cm catheter with 24cm fenestrated section, large diameter holes to maximise drainage and reduce occlusion risk.
- Complete with polyester cuff to promote rapid tissue in-growth, minimising infection risk and aiding catheter security.

SAFETY VALVE WITH UNIQUE LOCKING SYSTEM:

- Prevents leakage of air or fluid.
- Prevents accidental disconnection during drainage.
- Prevents accidental loss of valve caps, minimising contamination.

SIMPLE ORDERING:

- All the equipment needed for the drainage procedure and subsequent re-dressing is provided in a single pack: Product Code: **R54400 Rocket® IPC Dressing Pack & Bottle**
Further supplies can be obtained directly from **NHS Logistics** using the codes below.
- **The packs can be ordered on prescription tariff - Part IXA appliances: Indwelling Pleural Catheter Drainage Systems, Rocket IPC Drainage Pack & Bottle.**



Description	Code	NHS Code
Rocket® IPC Dressing pack & Bottle Available on Prescription Please prescribe drainage pack and bottle in Multiples of 5.	R54400	FET1544
Rocket® IPC Catheter Insertion Set	R54400-16-40	FET1545
Rocket® IPC Pre Evacuated Bottle Available on Prescription Please prescribe drainage pack and bottle in Multiples of 5.	R54410	FET1546
Rocket® IPC Drainage Line	R54410-00-DL	FET1547

