


Rocket[®] Fetal Blood Sampling Kits

Rocket[®] FBS Kits offer simple, safe and cost effective solutions for routine fetal blood sampling.

- **RELIABLE SAMPLING:** Developed in conjunction with the world's leading manufacturers of blood gas analysers to ensure compatible tubes and components.
- **MAXIMUM ILLUMINATION** with minimum of fuss. All Rocket[®] FBS Kits come complete with the revolutionary **Rocket[®] Amnilume™ Single Use Light Source**.
- **COMPREHENSIVE KITS:** each designed to complete up to 5 individual FBS samplings over the delivery period. Each kit is ready for immediate use to minimise delay when time is precious.
- **REDUCE RISK AND STRESS** caused by incompatible sampling tubes and missing or worn instruments frequently found in reusable trays.
- **EFFICIENT AND COST EFFECTIVE** when compared to reprocessing and replacing worn or lost reusable instruments.



Endoscope & Obturator
6 x 12" cotton swab sticks
5 x gauze swabs
1 x Disposable sponge holder
2 x Rocket[®] FBS Blades

1 x Rocket[®] FBS Blade Holder
1 x Yellow Soft Paraffin BP
1 x capillary tube holder
10 x capillary tube stoppers
Tray wrapped in waterproof drape

Packed individually in cartons of 10 units

	R57012	R57013	R57018	R57024	R57026	R57027
Amnilume™ Disposable Light Source	✓	✓	✓	✓	✓	✓
6 x swab sticks	✓	✓	✓	✓	✓	✓
5 x gauze swabs	✓	✓	✓	✓	✓	✓
Disposable Rampley Forceps	✓	✓	✓	✓	✓	✓
3 x long heparinised tubes, glass	Straight	Angled	N/A	N/A	N/A	N/A
5 high heparin plastic capillary tubes	75µl	75µl	55µl	85µl	100µl	170µl
3g Yellow soft Paraffin	✓	✓	✓	✓	✓	✓
Angled FBS Blades	2	2	2	2	2	2
Sampling Wand Version available	✓	✓	✓	✓	✓	N/A

