


Rocket® KCH™ Fetal Bladder Drain

The Rocket® KCH™ Catheter is intended for use in fetal urinary tract decompression following diagnosis of fetal obstructive uropathy in foetuses of 18-38 weeks gestation.



- **TWIN PIGTAIL** catheter with 6 holes for reliable drainage and metal tip markers for easier post insertion position monitoring.
- **MRI COMPATIBLE:** The catheter and marker rings are compatible with MR Imaging for simplified post insertion monitoring.
- **PTFE GUIDEWIRE** supplied to ease loading of pigtail catcher into introducer and minimises the risk of kinks
- **REUSABLE INTRODUCER SET** for secure, insertion and reliable placement.



Description	Code
Rocket® KCH™ Fetal Bladder Drain Set with 0.035" x 80cm PTFE guidewire, supplied sterile, pouch packed, each	R57405
Rocket® KCH™ Reusable Introducer Set: trocar, cannula, long & short pusher rods	R29100

