

# MityVac<sup>®</sup>

## Obstetric Vacuum Cup System

### R28134: Mityvac™ Hand Vacuum Pump:

- Trigger vacuum release gives the clinician complete vacuum control.
- Precision gauge is four colour-coded for easy interpretation and safety.
- Mini Grip handle provides maximum comfort for clinicians with smaller grip spans.
- Vacuum Range: 0-30 Hg, 0-76 cm Hg.
- Air Exhaust Rate: 36ml/stroke (at atmospheric pressure).
- Autoclavable: 135°C @ 2bar/3mins



### R57070-MUSH: M-Style™ Mushroom Cup with 4' Tubing & Filter

- Shaped to provide greater traction with gentle mechanical pull and excellent manipulation during vacuum assisted deliveries.
- Soft, contoured cup design allows for a stronger grip to minimise infant scalp trauma.
- Flexible stem bends into a 90° angle for ease of insertion whilst the Mushroom® Cup will maintain vacuum while bent at a 90° angle.
- Individually packed, sterile, for single patient use, 12 per carton.



### R57070-SELC: M-Select™ Mushroom Cup with 4' Tubing & Filter

- Advanced, versatile cup, providing optimum traction and manipulation during delivery.
- Smooth, low profile and string handle which allows for easier placement in Occiput Posterior (OP) and Occiput Transverse (OT) positions.
- Also suitable for low Occiput Anterior (OA) and Outlet presentations.
- Individually packed, sterile, for single patient use, 12 per carton.

