

Rocket® LLETZ Loops



The Rocket® LLETZ Loop Range was the world's first single use loop electrode designed specially for LLETZ. 20 years on the range is still in constant production. It continues Rocket Medical's tradition of providing exceptional value, reliability and quality

- **OVERCOME THE PROBLEMS** of high cost, cross infection risks and health and safety hazards commonly associated with reprocessing reusable LLETZ electrodes.
- **Rocket® LLETZ LOOPS** are designed to deliver higher cutting performance with less thermal damage to the sample.
- **2.4mm (3/32") x 10cm shaft** to fit all standard UK, US and European diathermy systems
- **DELTA™** style loops available for increased wire rigidity and increased control of penetration.

Rocket® LLETZ Range	Size		Code
Rocket® LLETZ Diathermy Loop, for single patient use, packed sterile in cartons 25 units	Small	14mm	R57800
	Medium	18mm	R57801
	Large	22mm	R57802
	Square	10mm	R57803
Rocket® DELTA™ LLETZ Diathermy Loop, for single patient use, packed sterile in cartons 25 units	Medium	18mm	R57851
	Large	22mm	R57852
Rocket® LLETZ Ball Electrodes, for single patient use, packed sterile in cartons 25 units	Standard	5mm	R57804
	Small	3mm	R57807
Rocket® LLETZ Diathermy Extension, for single patient use, individually packed sterile		50mm	R57818
Rocket® NETZ™ Needles, for single patient use, packed sterile in cartons 25 units		15mm	R57860

