

# Blunt Dissection Chest Drainage Insertion Pack

The **Rocket® Blunt Dissection Chest Drainage Insertion Pack** is available for those instances when it would be inappropriate or difficult to use a small bore seldinger chest drain such as a traumatic incident or the drainage of blood.

For further information please visit our website: [www.rocketmedical.com](http://www.rocketmedical.com) or contact our Customer Services Team to arrange a demonstration from one of our Clinical Sales Specialists.

## Blunt Tipped Catheter:

Used in conjunction with a stiffener so the catheter can be inserted in a more controlled manner.

## Curved Roberts Forceps:

To allow controlled dissection of the intercostal muscles over the rib.

## Radiopaque Stripe:

To allow accurate checking of catheter position.



## Insertion Pack:

Including a gown, drapes, suture and dressings.

**Simple Ordering** packs can be ordered direct from  
**Rocket Medical on 0191 419 6988**

Description	Size Fg	Code
Curved Roberts forceps, skin knife, guide wire, catheter with radiopaque stripe & stiffener, syringe, connectors, suture, ward procedure pack including gown, drapes, swabs, sponges, 10ml syringes, gallipot, wound dressing. Sterile, for single patient use in cartons of <b>5</b> units.	<b>26Fg</b>	<b>R54549-26-PK</b>
	<b>28Fg</b>	<b>R54549-28-PK</b>

