


Rocket® Water Trap Set

for use with Craft™ Suction Pumps

The Rocket® Water Trap Set for use with Craft™ Suction Pumps offers a simple, cost effective solution to contamination risks and reuse problems.

- **COMPLETE SET:** Provides the complete bottle and tubing set in single, sterile pack for one time patient use.
- **GLASS TRAP BOTTLE:** To maintain clear view in the event of potential pump contamination with aspirate. Designed to withstand >450mmHg vacuum.
- **REMOVES CONTAMINATION RISK** from irregular or uncontrolled decontamination of reusable bottle, bung and angle tube set.
- **PATIENT TUBE SET:** 2m patient tube set with integral luer lock connector to ensure secure, reliable connection to the aspiration needle set.
- **SUITABLE FOR USE:** Compatible with all versions of **Craft™ Suction and Duo Vac.**
Also compatible with the *genX and IVF On-line Pioneer Pro Pumps.*



Ordering Information

Rocket® Water Trap Set: suitable for all Craft™ Suction Pumps
Sterile, for single patient use.

10 per carton

R57685

