

Rocket® Double Lumen Oocyte Aspiration Needles

Rocket® Double Lumen Oocyte Aspiration Needles are the latest generation of oocyte aspiration needles with new features to improve performance.

- **HIGH HARVESTING & FLUSH RATES.** Rocket® Double Lumen Needles have the highest harvest and flush rates of any comparable needle system.
- **SX™ SPIRAL ECHO BEVEL DESIGN** with high resolution echogenic surface preparation to give smoother, easier tissue penetration with extreme ultrasound definition.
- **TWO NON KINK TUBING SET STYLES:**
Standard: two limb tube set with luer fittings into clear silicone bung
Oxford: single limb tube set into clear silicone bung with 45° cannula connection (shown).
- **SHAPED FINGER GRIP** with bevel marker and enhanced surfacing to significantly increase comfort and control
- **FLUSHING LIMB STOPCOCK** fitted as standard to prevent accidental emptying of flushing syringe on oocyte aspiration



Ordering Information	Size/Length	Tube Set	Code
Rocket® Double Lumen Oocyte Aspiration Needles etched echo band, polyurethane tube set with clear silicone stopper for Falcon® tubes, pouch packed, supplied sterile, 10 per carton, for single patient use	16g x 32cm x 55cm	Standard	R57608
	16g x 32cm x 70cm	Oxford	R57610
	16g x 35cm x 70cm	Oxford	R57611
	16g x 35cm x 90cm	Oxford	R57611-00-90
Rocket® SX™ DL Oocyte Aspiration Needle SX™ spiral echo band, polyurethane tube set with clear silicone stopper for Falcon® tubes, pouch packed, supplied sterile, 10 per carton, for single patient use	16g x 33cm x 70cm	Oxford	R57710-SX-70

