


Rocket® ESL™ & EDL™ IUI Catheters

Rocket® ESL™ IUI Catheters from Rocket Medical provide the ideal solution for routine intrauterine insemination, whereas Rocket® EDL™ Double Lumen IUI Catheters, are designed for insemination where cervical access is complicated by excessive retroversion, anteversion or cervical stenosis.



Specifications:

ESL™

Outer 131 x 2.60 x 2.23mm. – polypropylene
 Inner: 176 x 1.58 x 0.71mm.
 45mm active length – polythene

EDL™

Outer: 133 x 2.60 x 2.23mm. - polypropylene
 Inner: 173 x 2.10 x 0.55mm
 40mm active length – polythene.
 Inner Sylet: stainless steel.

- **ESL™ SINGLE LUMEN**, dual diameter catheter with 5 x 1cm markings, atraumatic distal tip and dual sides eyes for easy, gentle insemination.
- **EDL™ DOUBLE LUMEN** catheter with integral stilet wire, 5 x 1cm marking, atraumatic distal tip and dual sides eyes for insemination where a malleable catheter is required to gain access.
- **COMPETITIVELY PRICED** to provide excellent value in routine IUI.

	Length	Code
Rocket® ESL™ IUI Catheters. Sterile, for single patient use, peel pouch packed in cartons 25 units	18cm	R57621
Rocket® EDL™ IUI Catheters Sterile, for single patient use, peel pouch packed in cartons 25 units	18cm	R57622

