


# Rocket® Endometrial & Mucus Sampling Syringes

- **Rocket® ENDOMETRIAL SAMPLER** (Blue) for the detection of endometrial carcinoma, monitoring of the endometrium in HRT and endometrial dating to assist diagnosis in infertility.
- **Rocket® MUCUS SAMPLING SYRINGE** (Red) for post-coital sample collection, mucus viscosity evaluations, detection of mucosal contaminants for diagnosis in infertility. Also indicated for sampling for STD culturing and removal of mucus plug prior to insemination or embryo transfer.
- **WELL TOLERATED:** Only 3.1mm OD for easy, low trauma insertion. Depth markers on RS57.509 to minimize risk of perforation
- **SMOOTH PISTON ACTION** with high vacuum for rapid effective sampling, minimizing cross contamination of sample.

**Specifications:**

**Endometrial Sampler:** blue piston  
3.1OD x 240mm x 1.8ml volume. Polypropylene outer tube – side eye,

**Mucus Sampler :** red piston  
3.1OD x 240mm x 1.8ml volume. Polypropylene outer tube – end eye,



Ordering Information	Pack	Code
<b>Rocket® Endometrial Sampling Syringe.</b> Sterile, for single patient use	Sterile per 25	<b>R57509</b>
<b>Rocket® Mucus Sampling Syringe</b> for single patient use	Sterile per 25	<b>R57508</b>
	Clinically Clean per 50	<b>R57507</b>

