



Rocket® BLUE Disposable Instrument Range

The **Rocket® Blue** Range of disposable instruments are designed to provide users with good quality, sturdy all-metal alternative to part plastic disposable instrumentation.

All instruments are marked with a non-toxic blue coating which is rapidly becoming an Internationally recognised standard to indicate metal instrument disposability to users.

User recognition of the disposability of instruments ensures that they are disposed of at point of use, which eliminates the inconvenience to Hospital Sterile Supply Facilities of having to receive and identify instruments that have been mistakenly returned.

The **Rocket® Blue** Range instruments are supplied sterile, single wrapped for single patient use.

Description	Qty Per Box	Qty Per Carton	Cat No
Dressing Scissor Str Sh/Bl 13cm	20	240	R50150-S
Dressing Scissor Str Sh/Sh 13cm	20	240	R50175-S
Iris Scissor Str 11cm	20	240	R50151-S
Iris Scissor Cvd 11cm	20	240	R50174-S
Spencer Wells Artery Fcps Str 12.5cm	25	50	R50154-S
Spencer Wells Artery Fcps Str 20cm	25	50	R50152-S
Mosquito Artery Fcps Str 12.5cm	25	50	R50156-S
Mosquito Artery Fcps Cvd 12.5cm	25	50	R50155-S
Dissecting Fcps Plain, Serrated, B.E 12.5cm	20	240	R50159-S
Dissecting Fcps T.O.E Plain, Serrated 12.5cm	20	240	R50160-S
Gilles Dissecting Fcps 1 x 2 Teeth 15cm	25	50	R50161-S
McIndoe Dissecting Fcps 15cm	25	50	R50162-S
B.P. Scalpel Handle No.3	20	240	R50157-S
B.P. Scalpel Handle No.4	20	240	R50158-S
Mayo Hegar Needleholder 14cm	25	50	R50173-S
Probe With Eye 12.5cm	20	240	R50170-S
Jobson Horne Probe 18cm	40	80	R50169-S
Foerster Sponge Fcps 25cm	20	40	R50164-S
Lister Sinus Fcps 12.5cm	25	50	R50163-S

Call our Export Customer Services Team on +44 1923 651404

exportsales@rocketmedical.com